# Instructions for Use

## Instructions for Use / Handling Instructions for the Modular Osteosynthesis System

### IMPORTANT PRODUCT INFORMATION

Please read prior to clinical use!

Our products are manufactured according to high quality standards. The modular osteosynthesis system is made of titanium / titanium alloys and implant steel alloys complying with ASTM standards for surgical implant material.

#### DESCRIPTION

The modular osteosynthesis system consisting of plates, screws, instruments and cassettes is intended for the osteosynthetic correction of alveolar ridge defects. The surface of the implants is chemically passive and the material is antimagnetic. The implants are delivered non-sterile. No liability is accepted for products that have been modified in any way, used for a purpose other than specified or used or handled improperly.

#### NOTICES ON PROPER USE

Implants are precision-mechanical products and may only be used for their intended purpose by specially trained personnel. The proper selection and application of the products is the sole responsibility of the competent surgeon. Be aware that correct selection of the product is extremely important! Implantation at the correct anatomical location is important as well. Each implant should be selected according to the following criteria: bone defect to be corrected – state of health of the patient – proven tolerance of the implant material.

Overload or wrong selection can lead to screw or plate loosening or fracture and consequential loss of function. The surgeon is required to inform the patient on the advantages and disadvantages of the modular osteosynthesis system and the surgical intervention. In particular, point out that, due to the limited strength of the implants, overload must be avoided and that non-compliance may adversely affect the bone healing process.

It is extremely important to handle the products with care. This especially applies to transportation, cleaning, sterilization and storage. If plates need to be shaped during application, such adaptation must be limited to the absolute minimum required. Never bend the plates to and fro repeatedly! Excessive bending or correction can lead to cracks or subsequent implant fracture. Improper handling and care, as well as non-intended use, can lead to wear, failure or destruction of the implants. Moreover, implants can corrode or become functionally impaired if they are allowed to get in contact with aggressive substances.

All components of the modular osteosynthesis system must be from our own production and belong to the respective individual system. The components have been matched to each other during their development and production. The use of products from other manufacturers and/or systems can lead to unpredictable risks and therefore is not permitted.

Implants are intended for single-use only. Reuse is not permitted even if the products seem to be fit for reuse after application. Be aware that multiple use of an implant voids any warranty given by the manufacturer!

The lot number specified on our label provided on the packaging allows retracing the entire production process down to the raw material used. This lot number must be archived.

#### POTENTIAL COMPLICATIONS AND SIDE EFFECTS

Complications can arise in the following cases: 1) Patients who are unwilling or unable to follow instructions during the postoperative phase, e.g. due to their mental or neurological condition. 2) Patients with insufficient or low-quality bone tissue, circulatory disorders or latent infections. 3) Hypersensitivity to the implant materials used, i.e. allergic reaction to foreign bodies. In such cases, it is mandatory to perform appropriate tests prior to implantation. 4) Tendon ruptures caused by the implant (due to erosion). 5) Infections or inflammations. 6) Product failure due to fracture or deformation. 7) Bone deformation and refracture. 8) Poor fracture healing.

Activities involving high physical loads should be avoided after the operation.

#### STORAGE

The products of the modular osteosynthesis system must be stored in a dry and dust-protected place. The packaging must be removed prior to cleaning, disinfection and sterilization and the product placed into a suitable sterilization box (such as the cassettes offered by USTOMED INSTRUMENTE GmbH & Co. KG).

#### INSPECTION

Our final inspection does not release the customer from its obligation to check the goods upon receipt. Be sure to inspect the implants for fissures and cracks or other damage prior to use. Never use damaged implants! Some of our implants are provided with a colored anodic coating. Minor color deviations are possible and do not affect the quality of the implants.

#### CLEANING / DISINFECTION / STERILIZATION

The implants and instruments must be cleaned, disinfected and sterilized prior to use. This may only be done by trained personnel and in compliance with the regulations in force at the time. Most commonly used sterilization methods can be used for this purpose. We recommend using steam sterilization (fractionated vacuum – pulsing vacuum / Pre-Vac – method). Be sure to observe the instructions provided by the manufacturer of your sterilizer with regard to correct sterilization time, temperature, pressure setting and final drying time. Parameters for steam sterilization (fractionated steam sterilization method): Temperature: 133°C – 135°C (270–275°F); standard cycle: 15 minutes (min. 5 minutes); final drying time: 8 minutes.

USTOMED INSTRUMENTE GmbH & Co. KG assumes no responsibility for the proper sterilization of the products by the user.

Should you have any queries or need additional information, please do not hesitate to get in touch with us.

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