




Competence in Instruments
MADE IN GERMANY

CATALOGUE 2025


ustomed.de

 **USTOMED®**
INSTRUMENTE

READY FOR EXCELLENCE



 **USTOMED[®]**
INSTRUMENTE
READY FOR EXCELLENCE



 **USTOMED[®]**
INSTRUMENTE

READY FOR EXCELLENCE

02	<i>SCALPELS AND PERIOTOMES</i>	P.	7-9
03	<i>SCISSORS</i>	P.	10-16
04	<i>TC SCISSORS</i>	P.	17-18
05	<i>MICRO SCISSORS</i>	P.	19-20
08	<i>DISSECTING AND TISSUE FORCEPS</i>	P.	21-25
09	<i>MICRO AND ATRAUMATIC FORCEPS</i>	P.	26-29
10	<i>COLLEGE FORCEPS AND MEASURING INSTRUMENTS</i>	P.	30-35
13	<i>HEMOSTATIC FORCEPS</i>	P.	36-39
17	<i>AUXILLARY CLAMPS AND FORCEPS</i>	P.	40-43
19	<i>NEEDLE HOLDER</i>	P.	44-45
20	<i>TC NEEDLE HOLDER</i>	P.	46-55
24	<i>RETRACTORS</i>	P.	56-68
27	<i>ELEVATORS AND RASPATORIES</i>	P.	69-75
31	<i>BONE RONGEURS</i>	P.	76-78
42	<i>MEASURING DEVICES</i>	P.	79-80
46	<i>SINUS CURETTES</i>	P.	81-84

48	<i>MOUTH WEDGES, GAGS & PACKING INSTRUMENTS</i>	P. 85-90
67	<i>BONE PREP INSTRUMENTS AND TOOLS</i>	P. 91-104
68	<i>HARD AND SOFT TISSUE AUGMENTATION SCREWS AND PINS</i>	P. 105-113
73	<i>ORGANIZATION AND STORAGE</i>	P. 114-117
79	<i>CLEANING AND CARE</i>	P. 118
80	<i>EXTRACTING FORCEPS ENGLISH PATTERN</i>	P. 119-134
81	<i>EXTRACTING FORCEPS WITH ANATOMICALLY SHAPED HANDLE</i>	P. 135-141
83	<i>ROOT ELEVATORS AND LUXATOREN</i>	P. 142-154
84	<i>MOUTH MIRRORS AND SYRINGES</i>	P. 155-161
85	<i>HAND INSTRUMENTS PERIODONTICS</i>	P. 162-176
86	<i>HAND INSTRUMENTS</i>	P. 177-186
87	<i>PROSTHETICS</i>	P. 187-192
89	<i>PROSTHETICS</i>	P. 193-198
90	<i>STORAGE, ORGANISATION AND STERILIZATION</i>	P. 199-206
91	<i>INSTRUMENT SETS</i>	P. 206-221

02

SCALPELS AND PERIOTOMES



02

Standard Scalpels and Blades

Blades, sterile, 100pcs.



02-040-011
Size 11, steril,
100 pcs.



02-040-012
Size 12, steril,
100 pcs.



02-040-015
Size 15, steril,
100 pcs.



02-040-124
Size 12d, steril,
100 pcs.



02-040-153
Size 15C, steril,
100 pcs.



1/1



1/1



1/1



02-030-030

Scalpel Handle,
No. 3

02-030-070

Scalpel Handle,
No. 7

02-030-150 G

Scalpel Handle,
15 cm, round shaped
smooth

02-030-150

Scalpel Handle,
15 cm, round shaped

02-046-022

BLADE-BUTLER,
for removing scalpel
blades from handle

Micro Blades, sterile, 25 pcs.



02-044-063
Size 63, steril,
25 pcs.



02-044-067
Size 67, steril,
25 pcs.



02-044-069
Size 69, steril,
25 pcs.



1/1

02-038-140
Micro-Scalpel Handle,
anatomical shape, 14 cm



1/1

02-795-000POT3
Periotome, size 1, straight,
complete, peek handle



1/1

02-796-000POT3
Periotome, size 1, angled,
complete, peek handle

03

SCISSORS



REYNOLDS

**03-155-145**

14.5 cm, delicate, curved

JAMESON

**03-185-155**

15.5 cm, delicate, curved

KILNER

**03-449-120**

12 cm, delicate, curved

METZENBAUM

**03-071-145**


14.5 cm, curved

**03-081-145**

14.5 cm, delicate, curved

03

Suture & Gum Scissors



You can find the TC Version of this Scissor in Chapter 04!

IRIS



1/1



1/1



1/1

11.5 cm, straight

03-320-115

Standard

03-320-115 SC

Super-Cut

03-316-115

Standard, large ring

11.5 cm, curved

03-321-115

Standard

03-321-115 SC

Super-Cut

03-317-115

Standard, large ring

13 cm, X-tra large ring

03-327-130

curved

03-327-130 SC

curved, Super-Cut

03-326-130

straight

JOSEPH



JOSEPH

03-440-140

14 cm, straight

03-440-140 SC

14 cm, straight, Super-Cut



JOSEPH

03-441-140

14 cm, curved

03-441-140 SC

14 cm, curved, Super-Cut

KELLY



KELLY

03-630-160

16 cm, straight

03-630-18

18 cm, straight



KELLY

03-631-160 SC


16 cm, curved, Super-Cut

03-631-180 SC


18 cm, curved, Super-Cut

03

Suture & Gum Scissors



STANDARD:
Standard



SUPERCUT (SC):

Special Cut with tothing in one scissor blade in order to prevent sliding of gum and suture material

FOX



1/2

03-420-130

13 cm, straight



1/1

03-421-130

13 cm, curved,
plain



1/2

03-417-130 SC

13 cm, curved,
Super-Cut

LA GRANGE



1/1

03-401-105

10.5 cm
Standard

03-401-105 SC

10.5 cm, Super-Cut

LOCKLIN

Gum Scissors

03-635-160

16 cm, angled, serrated



SPENCER



1/1

Ligature Scissors

03-720-130

13 cm



1/1

03-325-11511.5 cm,
angled

1/1

03-351-11511.5 cm,
very fine tips

03

Special Scissors

SURGICAL SCISSORS

BEEBEE

LISTER



03-027-145

14.5 cm, curved

1/1



1/1

03-744-105

Wire-Cutting Scissors,
10.5cm, straight, sharp



1/1

03-781-145

Bandage Scissors, 14.5 cm

UNIVERSAL



1/2

pointy/blunt

03-022-145

14.5 cm, straight

03-022-175

17.5 cm, straight

03-745-105

Wire-Cutting Scissors,
10.5 cm, curved, sharp


03-805-180

Universal Bandage Scissors,
18cm

04

TC SCISSORS





TUNGSTEN CARBIDE (TC):

Special Hard-Metal Inlays increasing life time and sharpness, of the Instrument

FOX

JOSEPH



04-417-130

13 cm, curved,
serrated, Carbide



04-440-140

14 cm, straight,
Carbide



04-520-115

11.5 cm,
straight,
Carbide

04-521-115

11.5 cm, curved, Carbide

05

MICRO SCISSORS



I'MOTION

**05-803-180**very fine tips,
14.5 cm

curved, round

05-809-145

14.5cm

05-809-180


18cm

**05-621-145**very fine tips,
14,5 cm**05-801-145**

very fine tips, 14.5 cm

08

DISSECTING AND TISSUE FORCEPS





1/1



1/1



1/1



1/1

Dissecting forceps,
straight

08-020-145

14.5 cm

08-020-160

16 cm

Dissecting forceps,
straight, delicate

08-030-145

14.5 cm

08-030-160

16 cm

08-030-180

18 cm

Dissecting forceps,
curved, delicate

08-033-145

14.5 cm

08-033-160

16 cm

Dissecting forceps,
cross-serrated

08-036-145

14.5 cm


08-036-160

16 cm

08-036-180

18 cm

ADSON



1/1

08-060-122

12 cm, del. Mod.



1/1

Dissecting forceps

08-060-120

12 cm

08-060-150

15 cm

SEMKEN



1/1

08-062-150

15 cm, x-fine



1/1

08-080-155

15.5 cm, straight



1/1

10-130-18018 cm, small,
Carbide,
atraumatic

GERALD

SEMKEN



1/1

Tissue forceps, straight,
1x2 teeth, v. delic.

08-235-145

14.5 cm

08-235-160

16 cm

08-235-180

18 cm



1/1

08-281-155

15.5 cm, straight.,
1x2 teeth



1/1

08-301-175

17.5 cm, straight,
1x2 teeth

ADSON



Tissue forceps

08-261-120

12 cm

08-261-150

15 cm

MICRO-
ADSON

x-fine

08-262-120

12 cm


08-262-150

15 cm

ADSON-
BROWN**08-267-120**

12 cm, 7x7 teeth

09 MICRO AND ATTRAUMATIC FORCEPS



DE BAKEY

INFO
for tissues and sutures



1.5 mm

1/1



2.0 mm

1/1



1/1

Atraumatic

09-520-150

15 cm

09-520-180

18 cm

Atraumatic

09-526-150

15 cm

09-526-180

18 cm

09-545-150

15 cm, angled,
wide 2 mm



1/1

09-230-180

anatomical, 18 cm, straight,
diamanted suture plate



1/1

09-267-180

surgical, 18 cm, curved,
1x2 fine teeth



1/1

09-265-180

surgical, 18 cm, straight,
1x2 fine teeth

DE BAKEY



1/1

09-270-180

18 cm, straight,
atraumatic

I'MOTION




1/1

09-280-180

COOLEY, 18 cm,
straight, atraumatic

10 COLLEGE FORCEPS AND MEASURING INSTRUMENTS



LONDON
COLLEGE

LONDON COLLEGE modified

LONDON COLLEGE
modified

10-825-150

15 cm


10-825-160

16 cm

10-825-180

18 cm

FLAGG

**10-831-160**

16 cm, serrated jaws

MERIAM

**10-841-160**

16 cm, serrated jaws

10-845-155

15.5 cm, smooth handles



1/1

Splinter foreceps, 10.5 cm,
serr.

08-452-105

straight

08-453-105

curved



1/1

10-823-160

Maxgraft-Bone-Ring
forceps, 16 cm,
for Ø 6 mm a. Ø 7 mm Ring



1/1

10-871-160

16 cm, with eye tips

PONTE



1/1

10-877-160

Bone-Block holding
forceps, 16 cm, Micro



1/1

10-883-150

Nerve broach pliers,
15 cm, with lock

MILLER



1/1

10-947-150

forceps for articulating
paper, 15 cm

**10-941-000**

IWANSON, caliper, meas.,
range 0-10 mm - only for
comparing measures -



1/1

10-942-150

Ruler, 15 cm - only for
comparing measures -

**10-943-000**

measuring caliper - only for
comparing measures -

13

HEMOSTATIC FORCEPS



HALSTED-MOSQUITO



curved, anatomic

13-039-125

12.5 cm

13-039-140

14 cm



straight, anatomic

13-038-125

12.5 cm

13-038-140

14 cm

ROCHESTER-PEAN

MICRO-MOSQUITO



13-029-120
12 cm, curved,
anatomic



13-201-140
14 cm, curved,
anatomic

HALSTED-MOSQUITO

13-031-120

12 cm, curved, 1x2 teeth



1/1

OCHSNER-KOCHER



1/1

MICRO-MOSQUITO

HALSTED-MOSQUITO

13-043-125

12.5 cm, curved, 1x2 teeth

13-043-140

14 cm, curved, 1x2 teeth

13-042-125

12.5 cm, straight, 1x2 teeth



1/1

13-215-140

14 cm, curved, 1x2 teeth

17 AUXILLARY CLAMPS AND FORCEPS



GROSS MODIFIED



1/1

17-013-180

Sponge holding Forceps,
curved, with interlock

BACKHAUS



1/1

Towel Clamp, 8 cm

17-251-110

11 cm

17-251-080

8 cm



17-271-120

Towel Clamp,
paper/towels, 12 cm

JONES



17-227-090

Towel Clamp, 9 cm,
curved

STIEGLITZ



1/1

17-717-135

Splinter forceps, angled,
diamond-coated



1/1

Forceps, for seizing and
steril. of instruments, 20 cm

17-793-200

20 cm

17-790-145

14.5 cm

19

NEEDLE HOLDER





19-170-160

16 cm



19-300-140

14 cm


19-300-170

17 cm

20 TC NEEDLE HOLDER

TC = TUNGSTEN CARBIDE

An extremely durable metal, that is harder than most common suture-needles! This inlay increases the longevity of our needle holders dramatically. The inlays are able to change. The gold-plated handle indicates TC inserts.



HALSEY



1/1

20-120-130

13 cm



3/0 - 5/0



1/1

20-160-150

15 cm



4/0 - 5/0

CRILE-WOOD



20-164-150
15 cm

3/0 - 5/0



20-164-180
18 cm

3/0 - 5/0

MAYO-HEGAR



1/1

20-220-160

16 cm

 2/0 - 4/0

20-220-180

18 cm

 2/0 - 4/0

MICRO-RYDER



1/1

20-320-150

15 cm


 5/0 - 7/0

20-320-130

13 cm

 6/0 - 7/0

MATHIEU MODIFIED



20-452-170	17 cm	3/0 - 4/0
20-452-140	14 cm	3/0 - 4/0
20-452-200	18 cm	3/0 - 4/0

LICHTENBERG



20-471-170	17 cm, curved	3/0 - 4/0
------------	---------------	-----------

HÖSEL



1/1

20-480-200

20 cm

3/0 - 4/0

LANGENBECK



20-392-160

16 cm

3/0 - 5/0

20-392-180

18 cm

3/0 - 5/0

TOENNIS



20-398-180

18 cm

3/0 - 4/0

20

MICRO- NEEDLEHOLDER



BARRAQUER



20-710-140

14 cm, straight,
round handle



20-710-180

18 cm, straight,
round handle



20-711-140

14 cm, curved,
round handle



20-711-180

18 cm, curved,
round handle



I'MOTION



1/1



1/1

20-752-180

18 cm, straight

INFO

- extremely light
- plasma-coated tips
- easy to clean

20-753-180

18 cm, curved



5/0 - 8/0

24

RETRACTORS



GILLIES

KILNER



24-421-150
15 cm, 1 fork, sharp



24-425-001
18 cm, size 1, small



24-425-002
18 cm, size 2, large

JOSEPH

GUTHRIE



1/1

24-469-002

15 cm, 2 forks, sharp



1/1

24-573-002

16 cm, 2 forks, 2 mm



1/1

24-573-007

16 cm, 2 forks, 7 mm




24-021-001

16 cm, single fork, sharp



24-021-002

16 cm, double fork, sharp



24-023-001

16 cm, single fork, blunt

24

Cheek Retractors

MINI- LANGENBECK

CAWOOD- MINNESOTA

GIERL



24-153-000

Cheek Retractor

24-154-160

16 cm

24-157-155

15.5 cm



84-239-002

GLUCKMANN
PET-Shield holder
double-sided, incl. each
1 x Shield right and
1 x Shield left



INFO

Special tool for the Advanced Socket
Shield Technique: „PET“ by Dr. Howard
Gluckman.

Helps you to secure the gingiva, when
preparing the buccal root during the
procedure.



24-158-002

PET Shield, right



24-158-001

PET Shield, left

MIDDELDORPF

HILGER



1/1



1/1

24-159-001

Attachment only

24-159-000

complete with handle




1/1

24-161-002

22 cm, size 2, medium

LANGENBECK




Wound Retractor
21 cm

KOCHER

INFO

The flat front gives a wider contact field and therefore a much better grip.




Wound Retractor
22 cm

OBWEGESER

INFO

The flat front gives a wider contact area and therefore a much better grip. The groove allows a more atraumatic approach ... of the mouth



Wound Retractor
24 cm

KOCHER-
TROSSBACH

3/4

24-179-103size 3, 55 x 11 mm,
bended upwards


NENTWIG



3/4

24-191-01414 mm, serrated,
bended upwards

USTOMED



INFO

Angled, for
laying patients

INFO

Especially for
seating patients

INFO

Ergonomic Handle

24-189-180

18 cm, angled

24-187-190

19 cm, straight

LASTER

BRAANEMARK



24-186-150

80 mm wide, 150 mm long



24-203-105

Special retractor for
wisdom teeth removal

INFO

Cheek and lip spreaders self-expanding,
plastic material, autoclavable up to 134° C.



1/1

24-395-001

Cheek and
Lip Spreaders,
children



1/1

24-395-002

Cheek and
Lip Spreaders,
adults

24

27

ELEVATORS AND RASPATORIES



WILLIGER

6 mm



1/1

27-741-160

16 cm,
curved

7 mm

wide
Modell



1/1

27-745-160

16 cm, curved,
wide Mod.



1/1

27-743-130

13 cm,
curved

6mm



1/1

27-749-006

17,5 cm,
6 mm wide



7 mm
27-749-007

27-749-007

17,5 cm,
7 mm wide



1/1

9 mm
27-749-009

27-749-009

9 mm



1/1

11mm
27-749-011

27-749-011

11 mm

FREER

**27-451-180**18 cm, double-ended,
sharp/blunt**27-451-180POT3**18 cm, double-ended,
sharp/blunt,
peek handle

PRICHARD

**27-459-175**

17,5 cm, double-ended,

27-459-175POT317,5 cm, double-ended,
peek handle**27-717-175**17,5 cm, double-ended,,
curved, small/wide**27-717-175P0**17,5 cm, double-ended,
peek handle

PARTSCH

MOLT

GLICKMAN

**27-711-024**

Nr.24 G, double-ended

**67-711-002**


Nr.24 G, double-ended

**27-721-190**19 cm, double-ended,
curved

ZUCCHELLI

**INFO**

black coated tips,
allow better visuability
under gingiva



85-814-175P0T3




Mod. Z, peek handle

85-812-175P0T3

Papillae elevator, 17 cm,
peek handle

85-812-175

Papillae elevator, 17 cm,
peek handle

**85-816-170P0T3**Mod. B., 17 cm,
PEEK handle, black**85-818-001P0T3**Tunneling knife,
micro, curved,
PEEK handle, black**85-818-002P0T3**Tunneling knife,
curved,
PEEK handle, black**85-818-175P0T3**Tunneling knife,
curved,
PEEK handle, black



1/1

85-816-175

Mod.B, 17.5 cm



1/1

85-816-175P0T3

Mod.B, 17.5 cm,
PEEK handle



1/1

85-808-175P0T3

Mod.W, 18.5cm,
PEEK handle

31

BONE RONGEURS



MICRO



1/1

31-319-145

14 cm, 25° curved

FRIEDMAN



1/1

31-321-140

14 cm, curved

LUER



1/1

31-379-150

15 cm, medium

BEYER

FRIEDMANN-delicat



1/1

31-315-150

curved




1/1

31-421-180

18 cm, curved

42

MEASURING DEVICES



BOWMAN

INFO

Biegsam,
geeignet für
Thermo-
desinfektor



1/1

Antrum probe, ball-point

42-652-000
size 00/0



42-652-012
size 1/2



42-652-034
size 3/4



42-652-056
size 5/6



42-652-078
size 7/8




42-652-300
sz.000/00



42-652-400
0000/000



IMPLANTOLOGY CALIPER



1/1

42-850-020

Implantology caliper,
20 mm, straight

46

SINUS CURETTES



KHOURY SEMI-SHARP



1/1



1/1



1/1



1/1

46-695-004

Sinus Curette, Fig. 4,
180 mm,
semi-sharp

46-695-004POT3

Sinus Curette, Fig. 4,
180 mm, semi-sharp
PEEK handle

46-695-003

Sinus Curette, Fig. 3,
180 mm,
semi-sharp

46-695-003POT3

Sinus Curette, Fig. 3,
180 mm, semi-sharp
PEEK handle

HOTZ



46-695-001POT3

Sinus Curette, Fig. 1,
185 mm, semi-sharp

46-695-002POT3

Sinus Curette, Fig. 2,
190 mm, semi-sharp


46-695-001

Sinus Curette, Fig. 1,
185 mm, semi-sharp

46-695-002

Sinus Curette, Fig. 2,
190 mm, semi-sharp

NENTWIG SHARP

**46-699-001P0T3**

Sinus Curette, Fig. 1,
180 mm, sharp,
peek handle

46-699-002P0T3

Sinus Curette, Fig. 2,
190 mm, sharp,
peek handle

46-699-003P0T3


Sinus Curette, Fig. 3,
190 mm, sharp,
peek handle

46-699-004P0T3

Sinus Curette, Fig. 4,
190 mm, sharp,
peek handle


48 MOUTH WEDGES, GAGS & PACKING INSTRUMENTS





48-091-001

Mouth Wedge, Size S
37 x 15 x 22 mm (B x H x T)



48-091-002

Mouth Wedge, Size M
36 x 16 x 23 mm (B x H x T)



48-081-001

with steel core, Size M, 2 pcs.,
37 x 11 x 23 mm (W x H x D)



48-081-002

with steel core, Size L, 2 pcs.,
42 x 14 x 26 mm (W x H x D)

48-149-200

Silicone tube, 2 m, for
mouth gag, sterilizable

MOLT



1/1

48-141-140

14 cm, for adults

DENHART



1/1

48-145-120

12 cm, fine ratchet

48-149-200

Silicone Tube, 2m,
sterilisable, compatible
with 48-141-140,
48-145-120, 48-149-000

USTOMED

1/1

48-149-000

Set, complete, with two
tongue depressors for
children and adults



1/1

48-149-001

Scoop for 48-149-000,
single part, for children



1/1

48-149-002

Scoop for 48-149-000,
single part, for adults

HABERER

BRUENINGS

TOBOLD



1/1

48-260-200

20 cm, 17 mm x 25 mm



1/1

48-255-190

19 cm



1/1

48-271-165

16,5 cm

LUNIATSCHEK



1/1



1/1

48-339-002

Gauze Packer,
17,5 cm, breit

1/1



1/1

48-351-165

Retraction thread packer,
16,5 cm, gezahnt

1/1

INFO

Very thin and flexible working ends
slightly twisted in themselves, thus
making possible to operate more
simply and more carefully than with
the conventional method using the
Heidemann spatula.

48-339-001

Gauze Packer,
17,5 cm, schmal

48-349-165

Retraction thread packer,
16,5 cm, smooth

48-353-165

Retraction thread packer,
plate form

67/68 BONE PREP INSTRUMENTS AND TOOLS



DEGIDI



1/1

67-693-012P0T3

DEGIDI Plugger, size 1/2,
PEEK handle, for utilizing
and condensing of bone
material



1/1

67-693-034P0T3

DEGIDI Plugger, size 3/4,
PEEK handle, for utilizing
and condensing of bone
material

INFO

*Special instruments for
filling and condensing bone
material around an implant.*

RHODES

**67-545-367P0T3**

Back-Action Chisel,
Fig 36/37,
PEEK handle

85-802-305P0T3

Bone Scraper,
PEEK handle


67-691-003P0T3

Filler and Plugger
for biomaterials,
PEEK handle

67-691-008P0T3

SALAMA Filler
and Plugger,
PEEK handle,
for utilisation and
condensing of bone
material

MINI-LAMBOTTE



67-555-004 S	67-555-006 S	67-555-008 S	67-555-010 S
13 cm, 4 mm, curved, scaled	13 cm, 6 mm, curved, scaled	13 cm, 8 mm, curved, scaled	13 cm, 10 mm, curved, scaled
67-554-004S	67-554-006S	67-554-008S	67-554-010S
13 cm, 4 mm, straight, scaled	13 cm, 6 mm, straight, scaled	13 cm, 8 mm, straight, scaled	13 cm, 10 mm, straight, scaled

PARTSCH



1/1

67-711-002

Raspatory/Sharp
Spoon,
size 2, 15 cm


**67-678-170**

Bone Grinding Forceps,
grinding plates from
Titanium

NENTWIG


67-558-444

Set, straight, scaled
consisting of:




67-558-001

F.1, straight, scaled




67-558-002

F.2, straight, scaled



67-558-003

F.3, straight, scaled



67-558-004

F.4, straight, scaled


Localized spreading of the narrow ridge



NENTWIG


67-559-444

Set, baj., scaled
consisting of:




67-559-001

F.1, baj., scaled




67-559-002

F.2, baj., scaled



67-559-003

F.3, baj., scaled



67-559-004

F.4, baj., scaled

Condensing: preparation of the implant
site when bone structure is compromised



USTOMED BONE MILL



Set, complete

67-682-002

USTOMED Bone Mill, with 1 Cutting Cylinder Size S (ca. 2.0 mm)

67-682-003

USTOMED Bone Mill, with 1 Cutting Cylinder Size M (ca. 3.0 mm)

67-682-004

USTOMED Bone Mill, with 1 Cutting Cylinder Size L (ca. 4.0 mm)

INFO

Customer's Choice



67-682-012

Cutting Cylinder,
size S, (ca. 2.0 mm Chips) for
USTOMED bone flakes



67-682-013

Cutting Cylinder,
size M, (ca. 3.0 mm Chips) for
USTOMED bone flakes



67-682-014

Cutting Cylinder,
size L, (ca. 4.0 mm Chips) for
USTOMED bone flakes



67-682-020

Collector for
USTOMED bone flakes

INFO

exchangeable
teflon inlays

67-636-002

Teflon Discs, for
67-636-210, 2 pcs.

INFO

Completely
made from steel

INFO

Most common
mallet for use
with pins



1/1

67-660-150

16,5 cm, 150 g,
Ø 20 mm



1/1

67-636-210

18,5 cm, 210 gr,
Ø 25 mm

OMBREDANNE



1/1

67-670-720

26 cm, 720 gr,
Ø 40 mm



Photo Copyright Dr. Ponte



1/1

67-680-000

Bone crusher, round

LUCAS

LUCAS - serrated

**67-797-085**double-ended,
Fig. 85**67-797-085P0**double-ended,
Fig. 85, PEEK handle**67-797-086**double-ended,
Fig. 86**67-797-086P0**double-ended,
Fig. 86, PEEK handle**67-797-087**double-ended,
Fig. 87**67-797-087P0**double-ended,
Fig. 87, PEEK handle**67-797-186P0T3**double-ended,
Fig. 86, serrated,
PEEK handle

DDR, Lukas Salbrechter



18 cm, double-ended
67-799-000
size 0



67-799-010
size 1



67-799-020
size 2



67-799-030
size 3



67-799-040
size 4



1/1



67-799-010P0
size 1



67-799-020P0
size 2



67-799-030P0
size 3



67-799-040P0
size 4



1/1

18 cm, double-ended
67-799-000P0
size 0




1/1

INFO

Suitable for internal Sinus Lift.

Stopper is adjustable seamlessly.



Osteotom-Set (Universal Bego_RS-RSX_Set)

Bestehend aus:

Lagerungsrahmen und 5 Osteotomen
Durchmesser: 3,0/3,75/4,1/4,5/5,5 mm
Skalierung: 7/9/10/12/13/15 mm

67-851-020

Bone Condenser, bayonet, concave, Ø 2.0 mm

USTOMED TITAN-PIN-SET

68-780-018

Titan-Pin-Set, consisting of:



SET INCLUDES:



TITAN-PINS

68-620-003

Titan-Pins, 3 mm, 10 pcs.

INFO

Suitable mallet

67-660-150 (see page 98)



APPLICATOR

- Single-piece instruments to collect the pins from the holding tray therefore the twisting of the working tip during the operation is excluded.
- Very ergonomic and comfortable Grip for a safe operation of the Titan-Pins.

68-650-000

Titan-Pin-Applicator, 15,5cm



HOLDING TRAY

- safe and easy to operate with single hand principle
- Box can contain up to 15 Titan-Pins

68-720-008

Box for 15 Titan-Pins

68-780-016

Bone-Fixation-Set, Titan
 Ø 1.2/1.4/1.6 mm, consisting of:

**INHALT:**

68-720-002	Box for Micro-Screws, Ø 1.0/1.2/1.4/1.6 mm
68-740-002	Screwdriver, 11 cm
68-760-102	Blade, cross, self-retaining, for Ø 1.0/1.2/1.4/1.6 mm Screws, 60 mm
68-510-014	Pilot drill, Ø 1.0 mm, 14 mm
68-512-014	Pilot drill, Ø 1.2 mm, 14 mm
68-112-008	Micro-Screws, cross, Ø 1.2 mm, 8 mm, 10 pcs.
68-115-004	Membrane Screws Ø 1.5 x 4 mm, selfcutting, Titan

MATCHING ITEMS:**Accessories:**

68-720-002	Box for Micro-Screws for Micro-Screws, Ø 1.0/1.2/1.4/1.6 mm
68-740-002	Screwdriver, 11 cm
68-760-106	Blade, cross, self-retaining, short, for Contra-Angle, Length: 27 mm
68-760-102	Blade, cross, self-retaining, for all Ustomed Screws, Length: 60 mm

Pilot Burs:

68-510-014	Pilot Burs, Ø 1.0 mm, 14 mm
68-512-014	Pilot Burs, Ø 1.2 mm, 14 mm
68-514-014	Pilot Burs, Ø 1.4 mm, 14 mm

Micro Screws (10 pcs. per Pack):

68-112-006	Micro-Screws, cross, Ø 1.2 mm, 6 mm
68-112-008	Micro-Screws, cross, Ø 1.2 mm, 8 mm
68-112-011	Micro-Screws, cross, Ø 1.2 mm, 11 mm
68-112-013	Micro-Screws, cross, Ø 1.2 mm, 13 mm
68-114-006	Micro-Screws, cross, Ø 1.4 mm, 6 mm
68-114-008	Micro-Screws, cross, Ø 1.4 mm, 8 mm
68-114-011	Micro-Screws, cross, Ø 1.4 mm, 11 mm
68-114-013	Micro-Screws, cross, Ø 1.4 mm, 13 mm
68-115-004	Membran-Screws Ø 1.5 x 4 mm, self-cutting

68-780-010

Bone-Fixation-Set,
Implantat Steel

**CONTENTS:**

- 68-720-002 Box for Micro-Screws, Ø 1.0/1.2/1.4/1.6 mm
- 68-740-002 Screwdriver, 11 cm
- 68-760-102 Blade, cross, self-retaining, 60 mm length, for 68-740-002
- 68-508-014 Pilot drill, Ø 0.8 mm, 14 mm
- 68-110-011 S Micro-Screw, Steel, Ø 1.0 mm, 11 mm, Implant-Steel alloy

MATCHING ITEMS:**Accessories:**

- 68-720-002 Box for Micro-Screws, Ø 1.0/1.2/1.4/1.6 mm
- 68-740-002 Screwdriver, 11 cm
- 68-760-106 Blade, cross, self-retaining, 27 mm length, for Contra-Angle Handpiece
- 68-760-102 Standardklinge, self-retaining, for all Ustomed Screws, Length: 60 mm

Pilot Burs:

- 68-508-014 Pilot Burs, Ø 0.8 mm, 14 mm
- 68-510-014 Pilot Burs, Ø 1.0 mm, 14 mm

Micro-Screws (10 pcs. per Pack):

- 68-110-008 S Micro-Screw, Steel, Ø 1.0 mm, 8 mm, Implant-Steel alloy
- 68-110-011 S Micro-Screw, Steel, Ø 1.0 mm, 11 mm, Implant-Steel alloy
- 68-110-013 S Micro-Screw, Steel, Ø 1.0 mm, 13 mm, Implant-Steel alloy

BOX FOR MICRO SCREWS AND ACCESSORIES

- Sophisticated and functional Design
- Integrated ruler for the exact determination of screw and graft lengths (removable)
- Newly-developed Material by Aerospace Engineering - THERMOSTABLE
- All parts are suitable for machine processing in the thermos-disinfector as well as for all conventional methods, thus corresponding to the latest „RKI - Guidelines“
- Collection tray for all Ustomed screw sizes diameter (1.0 mm to 1.4 mm)




SCREW DRIVER

- perfect handle length – therefore suitable for all hand sizes
- The ergonomic handle allows a comfortable and fatigue-free operation
- Self-holding mechanism for collecting screws





Photo: Dr. Jan Kielhorn



INNOVATIVE GEOMETRY OF THE MICRO SCREWS

- Precisely milled cross-slit screw head
- all screw types (1.0 mm bis 1.4 mm) are to be collected with the same blade
- Fast and safe collection of the self-holding micro screws
- Extreme friction between screw head and blade, prevents wear-and-tear
- Much easier removal of debris from the cross-slit (in comparison to inner four-edge or Torx)




Photo: Dr. Stephan Beuer, M.Sc.



TITANIUM MICRO SCREWS

- diameter: 1.2 mm / 1.4 mm
- length: 6 mm / 8 mm / 11 mm / 13 mm



IMPLANT-STEEL ALLOY MICRO SCREWS

- diameter: 1.0 mm
- length: 8 mm / 11 mm / 13 mm
- perfectly suited for very small bone grafts and 3D-Augmentations
- special alloy material prevents Osseointegration and therefore a dramatically reduced danger of breaking the screws in comparison to titanium screws of similar lengths





Photo: Jan Kielhorn



Photo: Jan Kielhorn



Photo: Dr. Stephan Beuer, M.Sc.




Photo: Dr. Stephan Beuer, M.Sc.



Photo: Dr. Stephan Beuer, M.Sc.



Photo: Dr. Stephan Beuer, M.Sc.



T-FIX MEMBRANE SCREW

- conical threat
- lenght: 4 mm
- self-cutting/self-drilling, no pre-drilling
- perfect for applications in strong cortical bone
- for PTFE and titanium-reinforced membranes



Ø 1.5 mm, 10 pcs.



68-115-004

UMBRELLA-TECHNIQUE BY DR. DR. SCHLEE



IMPLANT-STEEL ALLOY

- diameter: 1.2 mm
- length: 8 mm / 11 mm / 13 mm
- screw head diameter: 5 mm
- no Osseointegration



68-780-014

SCHLEE, Umbrella-
Technique-Set Ø 1.2,
consisting of:

**INHALT:**

- 68-720-010 Box for Umbrella-Technique, for Ø 1,2 mm screws
- 68-740-002 Screwdriver, 11 cm
- 68-760-102 Blade, cross, self-retaining, for Ø 1.0/1.2/1.4/1.6 mm screws, 60 mm
- 68-510-014 Pilot drill, Ø 1.0 mm, 14 mm
- * 6x Umbrella screws, each 2x 8 mm, 11 mm and 13 mm

Umbrella Screw (packed individually)

- 68-112-208 S SCHLEE, Umbrella-Screw, cross, steel, Ø 1.2 mm, 8 mm
- 68-112-211 S SCHLEE, Umbrella-Screw, cross, steel, Ø 1.2 mm, 11 mm
- 68-112-213 S SCHLEE, Umbrella-Screw, cross, steel, Ø 1.2 mm, 13 mm

Accessories:

- 68-720-002 Box for Micro-Screws, Ø 1.0/1.2/1.4/1.6 mm
- 68-740-002 Screwdriver, 11cm
- 68-760-106 Blade, cross, self-retaining, Winkelst., for Ø 1.0/1.2/1.4/1,6 mm screws, 27 mm



73

ORGANISATION AND STORAGE





73-134-024

Bur black
12x FG
12x CA



73-100-050

Needle and Drill sterilizing box,
perforated, 50 x 25 x 4.5 mm



H: 8 mm

73-100-065

Needle and Drill sterilizing box,
perforated, 69 x 52.3 x 9,6 mm



73-200-003

Jar, 50 mm Ø,
130 mm



73-251-250

Kidney bowl, 250 x 140 x 40 mm

1/2

INFO*User's favorite***73-290-025**Metal-Bowl, 25 cc
Ø 60 mm**INFO***Heavy built,
for a safe stance***73-295-000**

Graft-Bowl

**73-292-050**Metal-Bowl, 50 cc,
graduated inside
Ø 55 mm**73-294-125**

Metal-Bowl, 125 cc

**73-294-050**

Metal-Bowl, 50 cc Ø 75 mm

79

CLEANING AND CARE

INFO

For maintenance of instruments prior to sterilization, for all instruments with joints, scissors, forceps, etc.



79-140-050

Passivation Bath,
corrosion-ex



79-130-040

Instruments oil spray, 500 ml,
biocompatible, sterilizable

80 EXTRACTING FORCEPS ENGLISH PATTERN



**80-001-000 HP**

Fig. 1, upper incisors
and canines, wide

**80-002-000 HP**

Fig. 2, upper incisors
and canines

**80-334-000 HP**

Fig. 34, upper incisors
and canines, deep-gripping



1/1

80-007-000 HP

Fig. 7, OK, upper premolars, HP



1/1

INFO

Reaches under
the gingival
margin

80-335-000 HP

Fig. 35, OK, upper premolars, HP

80

English Pattern, non-slip



INFO

„Peak to the Cheek“
(est. in KL 10.02.2025)

1/1

80-017-000 HP

Fig. 17, OK, upper molars, right



1/1

80-018-000 HP


Fig. 18, OK, upper molars, left



1/1

80-067-000 HP

Fig. 67, upper third molars



1/1

80-067-001 HP

Fig. 67A, upper third molars, narrow



1/1

80-051-000 HP

Fig. 51, upper roots



1/1

80-051-001 HP

Fig. 51A, upper roots



1/1

80-051-012 HP

Fig. 51L, upper roots

WITZEL



INFO

The 45° angle allows these forceps to be used for upper and lower roots.

1/1

80-355-000

Special Forceps, very fine, upper and lower roots



1/1

80-397-000 HP

Fig. 97, upper root fragments

INFO

for most delicate root
fragments

STIEGLITZ

1/1

80-555-000


Special Forceps,
very fine for roots



1/1

17-717-135

Splinter forceps, angled,
diamond-coated



1/1

80-597-000

Fig. 97, upper
root fragments

INFO

Reaches under the
gingival margin



1/1



1/1



1/1

80-013-000 HP

Fig. 13, lower incisors

80-013-001 HP

Fig. 13A, premolars,
narrow

80-336-000 HP

Fig. 36, lower incisors
and premolars,
deepgripping



80-022-000 HP

Fig. 22, lower molars




80-022-019

Fig. 22S, lower molars



80-086-001

Fig. 86A, (Cowhorn)
lower molars



English Pattern, non-slip

80



80-079-000 HP

Fig. 79, lower third molars, with tips



80-079-001 HP

Fig. 79A, lower third molars





80-300-000

Special Forceps, upper and lower wisdom teeth

80

English Pattern, non-slip



80-222-051 HP

Fig. 22, 1/2 left,
lower wisdom teeth

80-222-052

Fig. 22, 1/2 right,
lower wisdom teeth



80-346-012 HP

Fig. 46L,
lower roots

80-033-000 HP

Fig. 33,
lower roots

80-033-001 HP

Fig. 33A,
lower roots

80-033-012 HP

Fig. 33L,
lower roots

STIEGLITZ



1/1

17-717-135

Splinter forceps, angled



1/1

80-546-012Fig. 46L, lower roots,
Wurzelz

1/1


80-555-000Special Forceps,
very fine for roots



KLEIN


Extracting Forceps for Children

80




80-860-000

upper incisors,
Front




80-861-000

upper incisors,
premolars



80-865-000

upper incisors,
molars



80-863-000

upper incisors,
roots



80-867-000

lower incisors



80-869-000

lower incisors, molars

80-860-006

complete Set, consisting
of: 4 forceps upper jaw
and 2 forceps lower jar

80-000-010

BASIC EXTRACTION FORCEPS SET
consisting of:



80-002-000HP
upper Front




80-007-000HP
upper premolars




80-017-000HP
upper molars



80-018-000HP
upper molars



80-067-000HP
upper 3rd molars



80-051-000HP
upper roots



80-013-001HP
lower Front/
premolars



80-022-000HP
lower molars



80-079-001HP
lower 3rd molars



80-346-012HP
lower roots

81 EXTRACTING FORCEPS WITH ANATOMICALLY SHAPED HANDLE





1/1

81-002-000 HP

Fig. 2, upper incisors and canines, deep-gripping



1/1

INFO

Reaches under the gingival margin

81-334-000 HP

Fig. 34, upper incisors and canines, deep-gripping



81-007-000 HP

Fig. 7, upper premolars



INFO

Reaches under
the gingival
margin

81-335-000 HP

Fig. 35, upper premolars, deep-gripping



81-017-000 HP

Fig. 17, upper molars, right



81-018-000 HP

Fig. 18, upper molars, left



1/1

81-067-000 HP

Fig. 67, upper third molars, bayonet



1/1

81-051-000 HP


Fig. 51, upper roots, bayonet



1/1

81-051-001 HP

Fig. 51A, upper roots, narrow, bayonet



81-013-001 HP

Fig. 13A, lower incisors,
narrow



81-013-000 HP

Fig. 13, premolars




81-336-000 HP

Fig. 36, lower incisors and
premolars, deep-gripping

INFO

Reaches under
the gingival
margin



1/1

81-022-000 HP

Fig. 22, lower molars



1/1

81-086-001

Fig. 86A, (Cowhorn)
lower molars



1/1

81-079-000 HP

Fig. 79, lower third molars,
with tips

83 ROOT ELEVATORS AND LUXATOREN



BEIN



1/1



1/1



1/1



1/1

83-100-002

straight, 2 mm

83-100-003

straight, 3 mm

83-100-004

straight, 4 mm

83-100-005

straight, 5 mm

83

Root Elevators

BEIN



83-154-004

gezahnt, 4 mm

1/1



1/1

83-111-025

curved backwards,
2.5 mm



1/1

83-111-035

curved backwards
3.5 mm



1/1

83-111-045

curved backwards,
4.5 mm



1/1

83-117-035

modified, bayonet,
3.5 mm

USTO-LUX



1/1

83-120-003delicate, straight,
3 mm

1/1

83-120-004delicate, straight,
4 mm

1/1

83-120-005delicate, straight,
5 mm

USTO-LUX



1/1

83-121-003delicate curved,
3 mm

1/1

83-121-004delicate curved,
4 mm

1/1

83-121-005delicate curved,
5 mm

X-TRACTION - Pen Luxators



1/1



1/1



1/1



1/1

83-130-003

3 mm

83-130-004

4 mm

83-130-005

5 mm

83-130-104

BERNARD, Spear,
4 mm

FLOHR

CRYER



1/1

83-151-003

angled, Fig. 3



1/1

83-151-004

angled, Fig. 4



1/1

83-151-007

angled, Fig. 7



1/1

83-151-008

angled, Fig. 8

HEIDBRINK



1/1

83-231-009

angled, Fig. 9



1/1

83-230-008

straight, Fig. 8



1/1

83-231-010

angled, Fig. 10

FLOHR



1/1

83-427-001

angled, Fig. 1



1/1

83-427-002

angled, Fig. 2

CRYER



1/1

83-445-028

angled, Fig. 28



1/1

83-445-027

angled, Fig. 27


WINTER



83-741-006

Fig. 12 R, cross handle

BARRY




83-761-321

Fig. 321, cross handle



83-741-005

Fig. 12 L, cross handle



83-761-320

Fig. 320, cross handle

MOD. WIEN

**83-783-002**rotating elevator, 2 mm,
bayonet**83-783-003**rotating elevator, 3 mm,
bayonet**83-783-004**rotating elevator, 4 mm,
bayonet

HEIDBRINK

INFO

root fragments to be loosen
with HEIDBRINK, and
removed with STIEGLITZ

STIEGLITZ



1/1



1/1



1/1



1/1

83-871-001

Fig. 1

83-871-002

Fig. 2

83-871-003

Fig. 3

17-717-135

Splinter forceps, angled,
diamond-coated

83-871-001PO

Fig. 1,
double-ended
peek handle

83-871-002PO

Fig. 2,
double-ended
peek handle

83-871-003PO


Fig. 3,
double-ended
peek handle



1/1

83-904-000

size 0, straight,
hollow handle



1/1

83-905-002

size 2, angled,
hollow handle



1/1

83-905-001

size 1, angled,
hollow handle



1/1

83-899-002

size 2, angled



1/1

83-899-001

size 1, angled



1/1


83-898-000

size 0, straight

84


MOUTH MIRRORS AND SYRINGES



INFO*suitable for thermodesinfection***84-100-018**


ring handle, with aspiration hooks, stainless

1/1

INFO*suitable for thermodesinfection***84-060-018**

self-aspirating, with aspiration hooks, ring handle, stainless

1/1

INFO*NOT suitable for thermodesinfection, only manual disinfection***84-080-018**

with aspiration hooks, ring handle, compact

1/1

PAROJECT (RØNVIG)



INFO

suitable for thermodesinfection



84-149-000

Syringe with automated-pressure regulation
for intra-ligamental use

84-161-001

Spare cannula,
Ø 2.0 mm, Luer-Lock



1/1



1/1



1/1



1/1



1/1

84-218-000anatomical
shape**84-220-000**solid,
round shape**84-225-000P0**

peek handle

84-225-000P7**84-224-000**hexagonal
shape**86-080-000**Anatomical
USTOMED handle,
for probes and
mirrors



1/1

84-224-004

Combined mirror
size 4, 22 mm

INFO
12 pcs./pack

3
20 mm

4
22 mm

5
24 mm

plane		84-253-004	84-253-005
high-reflection (approx. 20%)		84-265-004	84-265-005
plane, Rhodium (no double pictures)	84-273-003	84-273-004	84-273-005

INFO
User's choice

INFO
In contrast to conventional mirrors, Rhodium
frontsurface mirrors have the reflecting metal
coating applied to the surface of the mirror plate
resulting in a **sharply** defined **refraction-free** image.

INFO

suitable for standard
mirror handles



1/1

84-285-003P0

plane, Rhodium,
Ø 3 mm,
peek handle black



1/1

84-285-005P0

plane, Rhodium,
Ø 5 mm,
peek handle black



11/1

84-285-062P0

plane, Rhodium,
6 x 3 mm,
peek handle black

**84-285-103**

3 mm,
insert only

**84-285-162**

6 x 3 mm,
insert only

**84-285-105**

5 mm,
insert only



1/2

84-295-170

Mirror for buccal site



1/2

84-292-210

Mirror for lingual site



84-296-300

Mirror Set (3 pcs.)



1/1

84-290-225

Mirror for occlusal site

85

HAND INSTRUMENTS PERIODONTICS



CUMINE



1/1

85-071-001

Fig. 1



1/1

85-071-002

Fig. 2



1/1

85-071-003

Fig. 3



1/1

85-197-000

RHEINWALD lifter

85-197-000PO

peek handle

WHITE

 1.6 mm 1.6 mm 1/1		 2.5 mm 2.5 mm 1/1	
 1.8 mm  2.1 mm  1/1  1/1		 1.8 mm  1.5 mm  1/1  1/1	
85-321-017 Fig. 17, double-ended, 1,8 + 1,8 mm	85-321-018 Fig. 18, double-ended, 2,1 + 2,1 mm	85-327-128 Fig. 127/128, double-ended, 1,8 + 1,8 mm	85-327-130 Fig. 129/130, double-ended, 1,5 + 1,5 mm
85-321-017P0 Fig. 17, double-ended, 1,8 + 1,8 mm, peek handle	85-321-018P0 Fig. 18, double-ended, 2,1 + 2,1 mm, peek handle	85-327-128P0 Fig. 127/128, double-ended, 1,8 + 1,8 mm, peek handle	85-327-130P0 Fig. 129/130, double-ended, 2,5 + 2,5 mm, peek handle
85-321-005 Fig. 5, double-ended, 1,6 + 1,6 mm		85-327-126 Fig. 125/126, double-ended, 2,5 + 2,5 mm	
85-321-005P0 Fig. 17, double-ended, 1,6 + 1,6 mm, peek handle		85-327-126P0 Fig. 127/128, double-ended, 2,5 + 2,5 mm, peek handle	

UNIVERSITY
NORTH
CAROLINA



1/1



1/1



1/1



1/1

85-691-015

CP15, 1-15mm,
colour marking

85-691-015P0

CP15, 1-15mm,
colour marking,
PEEK handle

85-691-115P0T3

CP15, 1-15mm,
colour marking,
double-ended with
Fig. 1 probe,
PEEK handle

85-695-012P0

CP12 - 3-6-9-12
peek handle black

85-695-012

CP12,
colour marking,
3-6-9-12

85-697-000

PSI/WHO,
colour marking

85-697-000P0

WHO/PSI
peek handle black

NABERS

**85-691-012**

Q2N, probe,
colour marking,
3-6-9-12

**85-699-012**

CP12, modified,
colour marking,
with tall tip,
0,5-3,5-6,5-12,5

**85-695-011**

CP11,
colour marking,
3-6-8-11

**85-693-010**

CP10,
colour marking

TAYLOR

INFO

User's choice

**85-745-023**M 23, double-
ended, universal**85-757-204**204 SD, double-
ended, universal**85-735-023**T2/T3, double-
ended, universal**85-735-023P5**T2/T3,
peek handle, grey**85-749-067**H6/H7, slight
angulation**85-741-023**CI 2/3,
double-ended**85-741-023P5**CI 2/3,
peek handle, grey**85-755-204P5**204 S, double-ended
PEEK handle, grey**85-755-204**

204 S, double-ended

GRACEY

**85-781-056POT3**

Curette, 5/6, PEEK handle

**85-781-078POT2**

Curette, 7/8, PEEK handle

**85-781-112POT4**


Curette, 11/12, PEEK handle

**85-781-134POT1**


Curette, 13/14, PEEK handle

85-779-056P7T3Curette, 5/6, PEEK handle
pink, DLC coated**85-779-078P7T2**Curette, 7/8, PEEK handle
pink, DLC coated**85-779-112P7T4**Curette, 11/12, PEEK handle
pink, DLC coated**85-779-134P7T1**Curette, 13/14, PEEK handle
pink, DLC coated


GRACEY (ONE SIDED SHARP)

**85-781-034**3/4, double-ended
Stahlgriff**85-781-056**

5/6, double-ended

85-781-056 POT35/6, double-ended
PEEK handle**85-781-078**


7/8, double-ended

85-781-078POT25/6, double-ended
PEEK handle**85-781-012**

1/2, double-ended

85-781-091

9/10, double-ended

**85-781-112**

11/12, double-ended

85-781-112 POT411/12, double-ended,
PEEK handle**85-781-134**

13/14, double-ended

85-781-134POT113/14, double-ended,
PEEK handle

COLUMBIA
UNIVERSAL

85-787-004

4L/4R, double-ended

M 23A
UNIVERSAL

85-791-023P9T0

Curette M23A, universal
peek handle orange

LANGER TITAN

**INFO***Titanium tips for
cleaning of implant
surfaces*

85-795-056P9T0

Curette, 5/6, double-ended,
sharp on both sides, blue
tips, orange PEEK handle

**85-719-067P5**Scaler, H6/H7,
peek handle grey**85-755-204P5**Scaler, 204 S,
peek handle grey

PERIO Set Gracey

**91-321-006**

PERIO Set Gracey:
 85-781-056P0T3, Fig. 5/6
 85-781-078P0T2, Fig. 7/8
 85-781-112P0T4, Fig. 11/12
 85-781-134P0T1, Fig. 13/14
 85-719-067P5, Scaler H6/H7
 85-755-204P5, Scaler 204S
 Inkl. washtray

**INFO**

With DLC-coating more durable sharpness than conventional curettes and scalers.

91-321-005

PA Set Gracey Advanced (DLC-coated)

Perio Set Gracey Advanced,
consisting of:

85-779-056P7T3,
85-770-078P7T2,
85-779-112P7T4,
85-779-134P7T1,
85-729-204S

ZUCCHELLI

INFO

*black coated tips, allow
better visuability under
gingiva*



1/1

85-814-175POT3

 Mod. W., peek handle


1/1


85-812-000POT3

17 cm, PEEK handle

1/1

85-812-175


Mod. W.

**85-816-175**

Mod. B, 17 cm

85-816-175P0T3Tunnelling knife,
micro, straight,
PEEK handle**85-808-175P0T3**PAR Dissector,
Mod. W.,
PEEK handle**85-808-175**PAR Dissector,
Mod. W., 17.5 cm

KASAJ

**85-816-170 POT3**

Tunnelling knife,
straight and Mod. B,
Raspa, PEEK handle

85-818-001POT3

Tunnelling knife,
micro, curved,
PEEK handle

85-818-002POT3

Tunnelling knife,
curved, PEEK handle

85-818-175 POT3

Tunnelling knife,
curved, PEEK handle

ORBAN

**85-831-012P0T3**

Gingivectomy knife,
19 cm, anterior,
PEEK handle

KIRKLAND

**85-825-156P0T3**

Gingivectomy knife, double-ended,
PEEK handle

86 HAND INSTRUMENTS






1/1

86-153-001

Size 1

86-153-001P0


Size 1, PEEK handle



1/1

86-153-017

Size 17



1/1

86-153-044


Size 44,
„Cow Horn“



1/1

86-157-001


Size 1, extra fine



1/1

86-157-017

Size 17, extra-fine



1/1

86-157-044

Size 44, extra fine,
„Cow Horn“

**INFO**

The „BLACK SPINE“ handle series features the special material PEEK (polyether ether ketone). This is a high-performance thermoplastic known for its exceptional mechanical, thermal, and chemical


properties. It is used in highly sensitive surgical procedures, such as spinal implants. This example has inspired us, for the name-giving of this extraordinary instrument series.

Why PEEK?

- + biocompatible and ideal for use in and on the human body
- + heat-resistant, acid-resistant, and X-ray neutral
- + non-degradeable, abrasion-resistant and corrosion-free

86-153-001P0

Dental Probe, Size 1,
single-ended,
PEEK handle



86-213-002

Furcation Probe,
unscaled
blunt tips

86-223-003

Fig. 3,
double-ended

86-211-002

double-ended,
extremely fine points



86-350-002


Fig. 2, double-ended



86-350-002P0

Fig. 2, double-ended,
PEEK handle

Plugger, double-ended, ball



INFO

Why PEEK?

- + biocompatible and ideal for use in and on the human body
- + heat-resistant, acid-resistant, and X-ray neutral
- + non-degradeable, abrasion-resistant and corrosion-free

HEIDEMANN

flexible/semi-sharp

Ø 2.0 mm
Stärke 0.55 mm



1/1

86-721-001P0

flexible, Fig. 1, PEEK handle

86-721-001

flexible, Fig. 1, round handle

Ø 2.5 mm
Stärke 0.55 mm



1/1

86-721-002P0

flexible, Fig. 2, PEEK handle

86-721-002

flexible, Fig. 2, round handle

Ø 3.0 mm
Stärke 0.55 mm



1/1

86-721-003P0

flexibel, Fig. 3, PEEK handle

86-721-003

flexible, Fig. 3, round handle

Spatulas (rigid/blunt), double-ended

Ø 1.8 mm
Stärke 0.4 mm



1/1

86-551-000P0Size 0,
PEEK handle**86-551-000**

Size 0

Ø 2.2 mm
Stärke 0.4 mm



1/1

86-551-001P0Size 1,
PEEK handle**86-551-001**

Size 1

Ø 2.0 mm
Stärke 0.5 mm



1/1

86-521-014

Size 14

Ø 2.5 mm
Stärke 0.5 mm




1/1


86-521-015

Size 15


Plugger (pear-shaped), double-ended



Plugger (plain), double-ended



Planisher, double-ended



2.0 mm
thickness: 0.4 mm

INFO

Tips are Titanium-
Nitrate coated



86-589-003


3TN

87

PROSTHETICS



R(H)EINSTOPFER



86-521-012P0T3

Size 2 retrograde
Plunger / Size 1
spatula,
PEEK handle

87-125-001P0T3

Fig 1, double-ended
PEEK handle

87-125-002P0T3

Fig 2, double-ended
PEEK handle

87-127-001P0T3

Fig 3, double-ended
PEEK handle

87-127-002P0T3

Fig 4, double-ended
PEEK handle

FAHNENSTOCK



1/1

Wax knife

87-290-130

13 cm

87-290-175

17,5 cm

AINSWORTH



1/1

87-420-000Rubberdam Punch,
17 cm




87-573-000

universal



87-570-000

straight

**87-591-000**

CROWN-BUTLER, semi-automated
Crown/Bridge remover

87-581-003

Crown Remover, complete
with 3 attachments

TOFFLEMIRE



87-344-001

1/1

TOFFLEMIRE,
Matrix Band, Fig. 1, 12 pcs.



87-344-002

1/1

TOFFLEMIRE,
Matrix Band, Fig. 2, 12 pcs.



87-344-003

1/1

TOFFLEMIRE,
Matrix Band, Fig. 3, 12 pcs.

TOFFLEMIRE-UNIVERSAL

TOFFLEMIRE-JUNIOR

TOFFLEMIRE-SENIOR



1/1

87-340-000

TOFFLEMIRE-UNIVERSAL,
Matrix Retainer



1/1

87-340-001

TOFFLEMIRE-JUNIOR,
Matrix Retainer



1/1

87-340-002

TOFFLEMIRE-SENIOR,
Matrix Retainer

89

PROSTHETICS



INFO

*For the atraumatic removal
of provisionals and bridges*

**89-120-140**

FURRER pliers, for
removal of provisionals
and copper ring

**89-115-140**

CROWN TRACTOR, complete
with 10pcs. Spare jaws and
10g adhesive powder
Manufacturer: Helmut Zepf
Medizintechnik GmbH

**89-115-020**

Exch. spare jaws,
synthetic,
Pack. 10 pcs.

**89-115-010**

Phiaal, cont. 10 g
adhesive powder



89-140-010

DEGIDI, CONO GRIP, silicone pads, small, 10 pcs.

DEGIDI



89-140-011

DEGIDI, CONO GRIP, upper jaw, delicate



89-140-023

DEGIDI, CONO GRIP, lower jaw



1/1

89-179-002

Diamond inserts,
only, 2 ea. pack.

**89-179-099**

Allan key,
only

**89-179-125**

BÖHM, telescope crown
pliers, diamond tips, incl.
Allen key and tips

ADERER

**89-554-140**

cross-cut jaws, 14 cm

**89-370-115**wire and clasp
bending pliers,
11.5 cm**89-260-160**UNIVERSAL-pin rough.,
cutt. and bend. plier,
16 cm



1/1

89-720-190

Plaster and laborat.spatula,
19 cm,
elastic bla.

90 STORAGE, ORGANISATION AND STERILIZATION



INFO

Storing Trays for an optimal Organisation of Hand Instruments during Cleaning, Sterilization and Storage.

Advantage: Tips of instruments are being saved and it speeds up the processing.

**90-340-004S1**

Storing Tray, Quarter-Size,
for 4 instruments, blue
177 x 56 x 27 mm

**90-340-005S1**

Storing Tray, Quarter-Size,
for 5 instruments, blue
179 x 86 x 26 mm

**INFO**

Fits the Demi Container
90-628-040 and 90-628-065

**90-340-008S1**

Storing Tray, Quarter-Size,
for 8 instruments, blue
178 x 128 x 26 mm

90-340-010S1

Storing Tray, Quarter-Size,
for 10 instruments, blue
178 x 128 x 26 mm

**90-340-020S1**

Storing Tray, Quarter-Size,
for 20 instruments, blue
178 x 250 x 26 mm

Tray item number +

Differently-coloured rubber inserts available. The colour definition is stated through the appendix of each item number:

S0 ● Black
S1 ● Blue
S2 ● Green
S3 ● Yellow
S4 ● Red

S5 ● Grey
S7 ● Violet
S9 ● Orange
S10 ○ White

90-626-040

Sterilization Container, Full-Size,
Pasteurian Loops, 4 cm

inner: 278 x 181 x 39 mm

outer: 313 x 190 x 54 mm

90-626-065

Sterilization Container, Full-Size,
Pasteurian Loops, 6.5 cm

inner: 278 x 181 x 54 mm

outer: 313 x 190 x 69 mm

90-626-085

Sterilization Container, Full-Size,
Pasteurian Loops, 8.5 cm

inner: 278 x 181 x 72 mm

outer: 313 x 190 x 86 mm

90-628-040

Sterilization Container, Demi-Size,
Pasteurian Loops, 4 cm

inner: 180 x 138 x 38 mm

outer: 211 x 145 x 56 mm

90-628-065

Sterilization Container, Demi-Size,
Pasteurian Loops, 6.5 cm

inner: 180 x 138 x 53 mm

outer: 214 x 145 x 68 mm

90-627-060

Sterilization Container, Slim-Size,
Pasteurian Loops, 4 cm

inner: 279 x 131 x 48 mm

outer: 310 x 139 x 66 mm

**INFO**

Sterilization Container
with permanent filter!

DENTAL STERILIZATION-CONTAINER SYSTEM WITH SAFETY-FILTER

Ustomed's answer to consistently increasing demands for risk minimization, environmental protection and handling comfort, is the newly-developed Easy-Stop System. The filterless mechanism stops microorganisms even more effectively than conventional filter systems and can be re-used as often as you like. There is no need for consumable purchasing. To open or close, one handling step is enough: Simply, press the lid down, and turn!

**THE ADVANTAGES SPEAK FOR THEMSELVES:**

The unique lid is almost indestructible:

- It is precision-milled from a solid block of aluminum. That way, the lid stays in shape and can withstand the daily office routines.
- Dirt- and germ-repellent
- Fast and residue-free drying during steam-sterilization
- High quality material for triple safety: adhere-free, pore-free and corrosion-free

INFO:

The two-piece aluminum insert in the lid serves as a germ barrier according to the principle of the „Pasteur Loops“. By forming a labyrinth of ring ribs, it therefore reliably separates microbes, germs and particles.

The meandering flow path achieves a decontamination rate of 99.997%



INFO

For storage of sterile
Sets

**90-620-040**

Sterilization Container, Full-Size, 4 cm,
one-filter system,
filter to be ordered separately

Innen: 281 x 184 x 39 mm
Außen: 327 x 193 x 54 mm

90-620-065

Sterilization Container, Full-Size, 6.5 cm,
one-filter system,
filter to be ordered separately

inner Dimensions: 281 x 184 x 56 mm
outer Dimensions: 327 x 193 x 71 mm

**90-618-100**

Sterilization Container, Full-Size, 10 cm,
one-filter system,
filter to be ordered separately

inner Dimensions: 275 x 180 x 92 mm
outer Dimensions: 315 x 190 x 105mm

90-618-130

Sterilization Container, Full-Size, 13 cm,
one-filter system,
filter to be ordered separately

inner Dimensions: 275 x 180 x 122 mm
outer Dimensions: 314 x 191 x 135 mm

INFO

Accessoires for
Container on page 205!



CUSTOMED TEFLON-FILTER

90-643-001

Teflon-Filter, for series
90-620-xxx

90-642-001

Teflon-Filter, for series
90-618-xxx



Ustomed Teflon-Filter

90-643-001

Teflon-Filter, for series 90-620-xxx

90-642-001

Teflon-Filter, for series 90-618-xxx

Marking Plates, including labeling

(max. 16 characters)



90-650-000

Red



90-650-001

Blue



90-650-002

Green



90-650-003

Yellow



90-650-004

Black



90-650-005

Silver



90-650-007

Purple



90-650-008

Orange



90-665-100

Sealing Plombs,
100 pcs.

or

90-665-101

Sealing Plombs,
1000 pcs.



90-654-000

Identification Label, plastic
Barcode, QR-Code Data-Matrix

or

90-656-000

Identification Label,
plastic Barcode 128



90-655-000

Plate for Stickers

90-655-000 for 1 pcs.

90-655-010 for 10 pcs.

90-650-025 for 25 pcs.



90-662-100

Indicator Labels,
100 pcs.

**90-716-001**

Wash Basket, Full-Size, MK4,
incl. 2 1/2 silicone inserts, cross-wise

inner Dimensions: 268 x 174 x 31 mm
outer Dimensions: 273 x 179 x 34 mm

90-716-002

Wash Basket, Full-Size, MK4,
incl. 2 silicone inserts, longitudinal

inner Dimensions: 268 x 174 x 31 mm
outer Dimensions: 273 x 179 x 34 mm

90-716-000

Wash Basket,
Full-Size, MK4, empty

inner Dimensions: 268 x 174 x 31 mm
outer Dimensions: 273 x 179 x 34 mm

**INFO**


suitable for sterilization
containers 90-626-xxx

90-716-050

Wash Basket, Demi-Size, MK4,
incl. 2 silicone inserts, suitable for
sterilization container 90-628-040

inner Dimensions: 172 x 133 x 31 mm

outer Dimensions: 177 x 138 x 35 mm




**90-716-025**


Wash Basket, Demi-Size longitudinal,
inkl. 2 silicone inserts, 2 baskets fit into 90-626-040

inner Dimensions: 274 x 88 x 35 mm

outer Dimensions: 268 x 83 x 31 mm

91 INSTRUMENT SETS





91-100-001

Essential Surgical Set

Constituent Parts:

02-030-030	Scalpel handle, No. 3
03-441-140	JOSEPH, suture-/gum scissors, 14 cm, curved
08-235-160	Tissue forceps, 16 cm, straight, 1x2 teeth, v. delic.
10-801-150	COLLEGE, tweezers, 15 cm, serrated jaws, standard
20-470-170	LICHTENBERG, needle holder, with TC, 17 cm
24-169-005	LANGENBECK, retractor, Fig. 5, 50 x 11 mm
27-717-175	Periosteal elevator, 17.5 cm, double-ended, curved
27-741-160	WILLIGER, periosteal elevator, 16 cm, curved
67-799-020	HEMINGWAY, bone curette, double, Fig. 2
84-224-000	Mouth mirror handle
86-153-001	Probe, Fig. 1, single-ended
90-716-001	USTOMED Wash basket, MK4 inkl. cross-wise-racks

Contents Essential Surgery Set 91-100-001



02-030-030



03-441-140



08-235-160



10-801-150



20-470-170



24-169-005



27-717-175



27-741-160



67-799-020




84-224-000



86-153-001



90-716-001



Constituent Parts:

91-101-001

Apical Surgery Kit

86-721-002	HEIDEMANN, spatula, flex., Fig. 2
86-521-002	Plugger, ball shaped, Fig. 2, double-ended
86-350-002	Cement spatula, Fig. 2, double-ended
84-285-062 P0	Root res. mirror, plane, rhodium, 6x3 mm
90-340-004 S1	Storing Tray, Quarter-Size, for 4 instruments, blue



DR. RICARDO KERN

We are proud to present the **Soft Tissue Tunneling Set**, developed by **Dr. Ricardo Kern**, one of Brazil's most renowned dentists. Dr. Kern specializes in the sophisticated design and reconstruction of soft tissue and has extensive experience in immediate implantation with provisional restorations, ranging from single-tooth cases to full-mouth rehabilitations.

This exclusive set includes five precisely crafted instruments specifically designed for tunneling soft tissue in dental procedures. With these high-quality tools, you can provide your patients with precise and comfortable treatment, whether for routine procedures or more complex surgeries. These instruments support you in making your work more efficient and effective.

Trust in the expertise of Dr. Ricardo Kern and expand your dental instrument collection with the **Soft Tissue Tunneling Set** to deliver outstanding results for your patients.



91-100-201

Set complete

Video Tutorial:



Constituent Parts:

85-818-101P9S0

KERN, Tunneling knife, Fig. 1, orange PEEK handle

85-818-104P9S0

KERN Papillae releaser „Rabbit Foot“/ Lingual Tunneling Knife, Fig. 4, orange PEEK handle

90-716-411

Wash tray „KERN“

85-818-102P9S0

KERN, Tunneling Knife, Fig. 2, orange PEEK handle

85-818-105P9S0

KERN Muscle releaser, Fig. 5, orange PEEK handle

85-818-103P9S0

KERN Papillae elevator, Fig.3, orange PEEK handle



Constituent Parts:

02-030-030	Scalpel handle, No. 3, flach
03-441-140	JOSEPH, suture-/gum scissors, 14 cm, curved
08-235-160	Tissue forceps, 16 cm, straight, 1x2 teeth, v. delic.
10-801-150	COLLEGE, tweezers, 15 cm, serrated jaws, standard
20-710-180	BARRAQUER, Micro-Needle Holder, TC, 18 cm, straight
24-169-004	LANGENBECK, retractor, Fig. 4, 40 x 11 mm
27-717-175	Periosteal elevator, 17.5 cm, double-ended, curved
27-741-160	WILLIGER, periosteal elevator, 16 cm, curved
67-799-020	HEMINGWAY, bone curette, double-ended, Fig. 2
73-100-065	Needle and Drill sterilizing box, perforated
84-224-000	Mouth mirror handle, hexagonal, s.s.
85-695-012	Periodontal probe, CP 12, colour marking, 3-6-9-12
86-153-001	Probe, Fig. 1, single-ended
90-716-001	USTOMED Wash basket, MK4 inkl. cross-wise-racks

91-100-002

Essential Implantology Set

Contents Essential Implantology Set



02-030-030



03-441-140



08-235-160



10-801-150



20-710-180



24-169-004



27-717-175



27-741-160



73-100-065



67-799-020



84-224-000



85-695-012



86-153-001



90-716-001



91-100-003

 Essential Micro Set

Constituent Parts:

02-038-140	Micro-Scalpel Handle, anatomical shape, 14 cm
05-809-180	Micro-spring handle scrs., 18 cm, curved, round handle
09-270-180	DE BAKEY, micro atraum. forceps, straight
20-710-180	BARRAQUER, TC Micro-Needle Holder, 18 cm, straight, round handle
85-691-015	UNIV. NORTH CAROLINA, period. probe, CP15, 1-2-3-4-5 skaliert
85-816-170 P0T3	Raspatory, Mod. B., 17 cm peek handle, Sperspitze
85-818-002 P0T3	Tunnelling knife, curved, peek handle
90-716-001	USTOMED Wash basket, MK4 inkl. cross-wise-racks



02-038-140



05-809-180



09-270-180



20-710-180



85-691-015



85-816-170P0T3



85-818-002P0T3



90-716-001



„I am using the I'Motion set for about a year. It is working perfectly with 6/0-7/0 threads. Thanks to it's small head and precise lock, I can work with it gently, furthermore, it is extremely rare to tear the thread, owing to it's rounded, smooth shape. It fits easily into narrow places.

PERFECT CHOICE NOT ONLY FOR PERIODONTAL PLASTIC MICRO- SURGERIES!

Dr. Dávid Botond Hangyási



I appreciate the ergonomic and lightweight design of these instruments providing a precise and

CONTROLLED HANDLING EVEN DURING LONG SURGERIES.

Prof. Dr. Adrian Kasaj, DDS, MSc, PhD

91-180-000

Mucogingival Plastic Surgery Set

05-803-180	I'Motion, Micro-scissors, 18 cm, curved
09-280-180	I'Motion, Micro-forceps, Cooley,
20-753-180	I'Motion, Micro-Needleholder, 18 cm, diam., curved alternatively 20-752-180
02-038-140	Micro-Scalpel Handle, 14 cm, anatomical shape
10-871-160	Suture forceps, 16 cm, with eye tips
85-691-115 P0T3	UNIVERSAL NORTH CAROLINA period. probe CP15+Probe
85-814-175 P0T3	Papilla elevator, Mod. B., double-ended
85-818-001 P0T3	KASAJ, tunnelling knife, micro, curved, peek handle black
85-818-002 P0T3	Tunnelling knife, curved, peek handle
85-818-175 P0T3	ORBAN Gingivectomy-knife, OR1/OR2
90-716-001	Wash Basket, Full-Size, MK4



**ULTRA-LIGHT
MICRO INSTRUMENTS**
THE NEW BENCHMARK IN SURGERY



02-038-140



05-803-180



09-280-180



10-871-160



20-753-180
or
20-752-180



85-691-115P0T3



85-814-175P0T3



85-818-001P0T3



85-818-002P0T3



85-818-175P0T3



90-716-001



DR. MARCO DEGIDI

Adjunct Professor, Postgraduate Program in Clinic Implantology and Biomaterials, University of Chieti. Adjunct Professor, Postgraduate Master in Clinic Implantology, University of Bologna.

Speaker in many national and international seminars, congresses and courses; many national and international publications.

Active member of AO, EAO, ICOI and SIO.

Private Practice in Bologna, restricted to implant surgery & prosthetic.



MR. GIANLUCA SIGHINOLFI

In 1989 Gianluca was awarded his diploma from the Villaggio del Fanciullo in Bologna, Italy. He immediately began work in an important dental lab in Bologna, later to become a partner.

He has had extensive specialist training and conducted courses both in Italy and internationally. Since 1996 Gianluca has collaborated with Dr. Marco Degidi and Prof. Piattelli on immediate loading research, contributing to the development of new surgery protocols and numerous articles on this subject that have been published in leading international reviews and journals.

Speaker in many national and international seminars, congresses and courses; many national and international publications.


surgery and prosthetics clinic, he knows exactly what he's talking about: perfection in practice means perfection for the patient.

TECHNIQUE

When used at implant level, these frameworks support extremely reliable temporary restorations, and rigidly stabilize immediately loaded implants, resulting in a dramatic improvement in implant success rates.


When used at abutment level with purpose-designed DENTSPLY Implants Manufacturing GmbH components, Intra Oral Welding makes it possible to manufacture extremely high quality and durable prosthetics with enormous precision.

The key to the Intra Oral Welding technique is an ingenious welding unit that was developed by the authors in collaboration with DENTSPLY Implants Manufacturing GmbH; the unit makes it possible to create a framework from a series of welds which join abutments to a connecting titanium bar. This quick and simple procedure is completely risk free for both surgeon and patient. There is absolutely no possibility of excessive heat and the procedure causes no discomfort of any sort to the patient. The very high quality of the weld obtained using this form of electrical resistance welding is comparable and often superior to traditional methods (laser, TIG, etc).



91-110-000

Degidi, WeldOne, IOW-Set



detailed information:

www.dentsplyimplants.com



Four implants are placed in the intraforaminal region.



Welding abutments in position.




Intra-orally welded framework.



Welded framework removed and reinforced.



Final restoration.



Post-op X-ray.

**Sharp Curver**

bends the bar into more acute angles for single abutments.

**Soft Curver**

bends the titanium bar along the curve of the arch.





Holder
grasps and holds the titanium bar firmly
in place during the welding process.




Cutter with tungsten carbide inserts
easily and efficiently cuts off excess
material.







 **USTOMED®**
INSTRUMENTE
READY FOR EXCELLENCE



 **USTOMED®**
INSTRUMENTE

USTOMED INSTRUMENTE GmbH & Co. KG

Bischof-Sproll-Straße 2
78532 Tuttlingen / Germany

Phone: +49 (0) 74 61 / 9 65 85 0

Fax: +49 (0) 74 61 / 9 65 85 65

info@ustomed.de

www.ustomed.de